

Warnings and Precautions:

This product utilizes a lateral flow chromatographic double-antibody sandwich method to detect labeled reporter molecules for Cas12/13 assays. Users simply need to label one end of the synthetic reporter molecule with Biotin and the other end with FITC (Fluorescein Isothiocyanate) or 6-FAM (6-Carboxyfluorescein).

This strip is designed to detect Cas enzyme cleavage products from nucleic acid amplification targets (such as PCR, LAMP, RPA/RAA).

Intended Use:

Detection of Cas enzyme cleavage products.

Packaging & Storage:

- **Packaging:** 50 strips/tube, moisture-proof aluminum foil bag.
- **Cat. No.:** JY0301
- **Storage:** Store in a dark, dry place at **4–30°C**.
- **Shelf Life:** **12 months**. After opening, reseal the tube cap tightly and store in a dry, dark place at 4–30°C (valid for 6 months)




Fig 1. Schematic Diagram of Disposable Nucleic Acid Test Strip

Sample Pad (Conjugate Pad): The area where the sample solution binds and initiates the detection process.

Test Area (Nitrocellulose Membrane): The region for observing and interpreting test results (containing Test and Control lines).

Absorbent Pad: Absorbs excess liquid during the assay to maintain flow.

Protocol:

1. Take out the required number of strips and mark the absorbent pad (Fig. 1). Each strip is for single use only.

2. Sample Preparation:

olf the cleavage product volume is 50–100 μL , it can be detected directly in a 200 μL PCR tube.

olf the volume is less than 50 μL , add ultrapure water to the PCR tube to bring the volume up to 50 μL and mix well.

3.Dilution: After CRISPR digestion of PCR, LAMP, RPA, or RAA products, dilute the cleavage product 2–20 fold (adjust based on product type and reaction system).

4.Detection: Insert the Sample Pad end (arrow end) of the strip into the reaction tube containing the diluted cleavage product (>50 μL). Ensure the liquid level does not exceed the top of the sample pad (Fig. 1).

5.Reading: Wait until the reading window is fully wetted (approx. 1–2 min; may take longer in low temperatures). Once the Control Line (C-line) appears, remove the strip, lay it flat, and read the result immediately.

6.Timing: Read results within 10 minutes after the C-line appears. Results read after 10 minutes are invalid.

7.Disposal: Record results and dispose of the strip safely in a sealed bag.

Result Interpretation:




Fig2. Cas12/13 Result Interpretation Guide for Cas12/13 Dedicated Nucleic Acid Test Strips

1. **Positive (+) :**

- Both Control Line (C) and Test Line (T) appear red.
- Only Test Line (T) appears red (Control Line C does not appear).
- *Note: Both scenarios indicate effective Cas12/13 cleavage activating the reporter.*

2. **Negative (-) :**

Only Control Line (C) appears red; Test Line (T) is not visible. This indicates Cas12/13 failed to cleave the reporter.

3. Invalid:

No bands appear on either C or T lines. This suggests the strip or reagents may be damaged/expired, or an operational error occurred. Retest with a new strip

Precautions & Safety:

- **Reporter Requirement:** Must be used with Biotin and FITC/FAM labeled reporters. Impure reporters (containing free Biotin or FITC) may cause false positives (red T-line in negative controls).
- **Negative Control Optimization:** Recommended reporter concentration range: **10–100 nM**. Perform a gradient dilution to find the concentration where only the C-line appears (no T-line) as the negative control baseline. Use this to calculate the optimal working concentration for your CRISPR system.
- **For Research Use Only (RUO).** Follow instructions strictly.
- **Storage:** Improper storage or use of expired products may lead to erroneous results. Use strips as soon as possible after opening to avoid moisture damage.
- **Single Use:** Do not reuse strips. Dispose of properly after use.