

Principle

RPApex Basic uses enhanced Recombinase Aided Amplification (RAA) technology for isothermal amplification at 39–42°C.

- Recombinase-primer complex invades target DNA template
 - Single-strand binding proteins maintain open-chain state
 - DNA polymerase initiates strand synthesis at 3' end
- **Exponential amplification achieved in 20–40 minutes**
- Ready-to-use lyophilized pellets for easy transport and storage

Primer Design Guidelines

- Length: 28–35 nt (longer than PCR primers)
 - Too short (<28 nt): reduced recombinase binding efficiency
 - Too long (>35 nt): secondary structures (hairpins, dimers) interfere with specificity
- **Concentration: optimize via gradient experiments (0.2–0.6 μM)**
- Amplicon length: 80–500 bp recommended
 - • GC content: 40–60% preferred
 - • Verify specificity via BLAST to avoid cross-reactivity

Intended Use

1. • For research use only (RUO)

1) • Amplification of DNA and RNA templates

- Applications: CRISPR detection, agarose gel electrophoresis analysis

Materials Supplied

- 48 reactions per kit
- Lyophilized pellets containing: recombinase, DNA polymerase, SSB, ATP, dNTPs, buffer components

2) Note: Mg²⁺ not included by default. Recommended final concentration: 14 mM/Test. Specify when ordering if required.

Storage and Stability

2. • **Unopened: Store at 25 – 30°C, dry and light-protected. Shelf life: 12 months (≥95% activity retention)**

- After opening: Tighten cap, continue storage at 25–30°C, dry and light-protected. Shelf life: 2 months

Component	Volume
Primer F (10 μ M)	1 μ L
Primer R (10 μ M)	1 μ L
Enzyme Lyophilized Pellet	1 pellet
Mg ²⁺	1 Mg ²⁺ pellet OR 1.25 μ L of 280 mM Mg ²⁺
Nucleic Acid Solution or Crude Lysate	2–4 μ L
Nuclease-free Water	Bring to 25 μ L total volume

3. Sample Preparation

Animal Samples (muscle or organ: spleen, tonsils, lymph nodes, kidneys, bone marrow):

- Cut 0.05 g or 3 mm³ tissue pieces into centrifuge tube
- Option 1: Add 200 μ L of 5 \times Lightning Reagent (BT0066), heat at 95°C for 5 min, use supernatant
- Option 2: Add 50–200 μ L of CoolFlash (BT0068) or StarFlash (BT0069) Reagent

4. Plant Samples (leaf, root, stem, or bud):

- Cut 2–3 mm² pieces into centrifuge tube

5. Add 50–200 μ L of CoolFlash (BT0068) or StarFlash (BT0069) Reagent

- Grind briefly, vortex to mix, use supernatant

6. Reaction Setup

Pipetting order: Negative control → Test samples → Positive control

7. Cap immediately after each addition to avoid aerosol contamination.

Mixing

8. Vortex at low speed or pulse mode for 3 – 5 seconds, OR

- Cap tightly and invert 3–5 times

• Centrifuge briefly and proceed to reaction

Incubation

• Place tubes at 39–42°C for 20–40 minutes

1. Post-Reaction Processing

1. Remove tubes from incubator immediately after reaction.
2. Add 50 μ L phenol/chloroform (1:1 v/v) to each tube.
3. Vortex vigorously for 10–15 seconds.
4. Centrifuge at 12,000 rpm for 1 minute to separate phases.
5. Pipette 10 μ L of supernatant (aqueous phase) for agarose gel electrophoresis.
6. Recommended gel concentration: 1.5%–2.0% (adjust based on amplicon size).

7. Visualize using appropriate staining (e.g., EB or GelRed).

Detection Sensitivity

